# **TOPO-REF-KERATOMETER** RT-7000

TOPOGRAPHY, AUTO REFRACTION & KERATOMETRY



- 4 systems 1 instrument
- Auto alignment + auto shot
- Colour touch screen
- Different topography maps
- Pupil + cornea ø measurementTear stability analysis system (TSAS)



# THE TOMEY TOPO-REF-KERATOMETER RT-7000



# QUALITY IN DETAIL

A unique combination of Topography, Autorefraction, Keratometry and TSAS (tear stability analysis system) in one - multi diagnostic replaces four devices with all their functions. The compactness of this instrument is its strength. It is therefore a perfect space and cost saving solution for you. Highly accurate easy handling makes working with the RT-7000 professional and quick.

#### **Colour touch screen**

The 6.4 inch coloured touch screen is used as operating monitor as well as for displaying all measured values. You can even move the unit in all directions by simply touching the screen. All commands can be done via touch screen.

#### Auto alignment + auto shot

The handling of the RT-7000 is very easy - it does almost everything by itself. Alignment and measurement are done automatically. You just roughly align the system towards the patient eye and the rest is taken care of by the instrument.

### Tear stability analysis system

Tear stability analysis system for analysing the tear film stability by using the light cone system in the RT-7000. That offers you several measurement and analysing functions to detect patients with dry eye indication.

#### Pupil + cornea ø measurement

cornea measurement bars to measure the individual diameter.

### Topography indices KRI + KAI

For immediate understanding of the cornea topographic structure we have implemented the topographic indices KAI (Kerato-Asymetry Index) and KRI (Kerato-Regularity Index). These values are highlighted in colour (green = normal, yellow = suspect, red = abnormal) to provide you a quick information about the corneal structure behaviour.

## **SPECIFICATIONS**

#### REFRACTIVE POWER MEASUREMENT

#### Spherical refractive power (S)

Measurement range -25.00 D to +22.00 D (at VD = 12.0 mm)

0.01 D, 0.12 D, 0.25 D Display unit

#### Cylindrical refractive power (C)

 $0 D to \pm 10.00 D (at VD = 12.0 mm)$ Measurement range

Display unit 0.01 D. 0.12 D. 0.25 D

#### Astigmatism axial (A)

Measurement range 0° to 180° Display unit

#### CORNEAL CURVATURE MEASUREMENT (K1, K2, AVG)

Measurement range 5.00 mm to 11.00 mm / 30.68 D

to 67.50 D (n=1.3375)

Display unit 0.01 mm

#### CORNEAL ASTIGMATISM AND AXIS (C, A)

Measurement range (C) 0 D to 10 D (n=1.3375)

Measurement range (A)

Measurement area cornea Ø 3.0 mm (at 8.00 mm corneal curvature)

#### **CORNEAL SHAPE MEASUREMENT**

Measurement area (at 8.00 mm corneal curvature) Normal measurement mode Ø 1.0 mm to 8.0 mm Special measurement mode Ø 0.9 mm to 7.0 mm

9 D to 100 D Display range

50 mm to 86 mm PD range Minimum pupil diameter Ø 2.2 mm Vertex distance 0 mm - 16.0 mm

#### **MAIN UNIT**

**Built-in printer** Thermal printer

Output External printer / LAN / USB


6.4" colour LCD Display Chin rest electr. controlled

#### **DIMENSIONS & ELECTRIC REQUIREMENTS**


**Dimensions WDH** 307 x 490 x 466 mm Approx. 20.0 kg Weight Voltage 100 VAC to 240 VAC Frequency 50/60 Hz

**Power consumption** 120 VA to 150 VA

# **DIMENSIONS**




PUPTIL




Pupil & cornea Ø measurement

Auto aligment

auto shot



Dual map



Post lasik normalized map

2012/04 - subject to change without notice

**TOMEY EUROPE TOMEY GmbH** Am Weichselgarten 19a 91058 Erlangen, Germany

Phone +49 9131 77710, Fax +49 9131 777120 Email info@tomey.de

**TOMEY ASIA-PACIFIC TOMEY CORPORATION JAPAN** 2-11-33 Noritakeshinmachi Nishi-ku, Nagoya 451-0051, Japan

Phone +81 52 581 5327, Fax +81 52 561 4735 Email intl@tomey.co.jp

