

LOGIC 65LG

Hi-tech ophthalmic equipment

Two instruments unit manual positioning of the table

Right and left solution

Motorized elevation 200 mm

Housing for handles instrument

*Wide range of accessories:
Chest of drawers with trial lenses, vertical or
horizontal counterbalanced arm, supports
for projector, video, ecc*

*Suitable for any model of chair of
our production*

*Suitable for easy access by
wheelchair users and for
sight visit to patients with
disabilities*



LOGIC 65LG

Hi-tech ophthalmic equipment

CONTROL PANEL



TRIAL LENSES DRAWER



CHAIR RAISING PUSH BUTTONS



MINIMUM INSTRUMENTS HEIGHT



Using the omnidirectional base BCRS applied to the chair, is possible to have easy access with a wheelchair to visit disabled patients

The logic unit, has a minimum height of the instruments table of 700mm with a raising of 200mm, allowing a better range of possibility for the sight examinations



Power



Projector



Auxiliary outlet



LED Lamp



LED Spot



Room Light



Up and down instruments table



Chair Raising, lowering and reset



Refractor arm opening/closing

LOGIC 65LG

Logic unit with wheelchair and 88MC chair on movable base BCRS.

TECHNICAL DATA

Input tension:	220 V 50/60 Hz 110 V 50/60 Hz
Power:	1060 VA
Led lamp:	12 V
Table heights:	min 700 mm max 900 mm
Weight only unit:	120 Kg
Counterweight:	45 Kg
Width:	~1600 mm
Depth:	1650 mm
Height:	1940 mm



POLYVISUS®

