




AJL BBG



The future of vision

AJL BBG Ophthalmic Solution

<i>Reference</i>	AJL BBG
<i>Model</i>	BBG ophthalmic solution at 0.025% W/V.
<i>Brand name</i>	AJL BBG ophthalmic solution.
<i>Definition</i>	Brilliant blue G ophthalmic solution, AJL BBG 0.025% W/V.
<i>Description</i>	AJL BBG 0.025% is a sterile, non-pyrogenic, isotonic solution, ready for use in vitreoretinal surgery. The BBG brilliant blue solution is capable of effectively staining the inner limiting membrane (ILM).
<i>Composition</i>	Trypan brilliant blue G sterile ophthalmic solution at 0.025% W/V.
<i>Indications / Properties</i>	<p><u>Indications:</u></p> <p>Greater affinity with the inner limiting membrane (ILM) and lesser with the epiretinal membrane, allowing an improved view of the ILM and facilitating the peeling in surgical operations for macular diseases, such as Macular Hole or Macular Fibrosis.</p> <p><u>Properties:</u></p> <ul style="list-style-type: none"> • Ready-to-use sterile solution: does not require mixing or filtering. • Excellent biocompatibility. • No adverse side-effects. • No side reaction such as photo induced cross linking of collagen fibres. • Easy removal of residual product from the eye.
<i>Sterilization method</i>	Steam.
<i>Supply</i>	Sterile, single unit. Contains a 1ml prefilled glass syringe and Luer lock 27G cannula.
<i>Shelf-life</i>	3 years

CE 0434

REV: 03