



AJL BLUE



The future of vision

AJL BLUE Ophthalmic Solution

<i>Reference</i>	AJL BLUE
<i>Model</i>	Trypan blue ophthalmic solution at 0.06% W/V.
<i>Brand name</i>	AJL BLUE ophthalmic solution.
<i>Definition</i>	AJL BLUE trypan blue ophthalmic solution at 0.06% W/V.
<i>Description</i>	AJL BLUE 0.06% is a sterile trypan blue solution, intraocular and biocompatible used in ophthalmic surgery of the anterior segment. This ophthalmic dye is used as a coadjuvant during capsulorhexis and minimises the risk of laceration in eyes with mature cataracts, with poor reflexes at the back of the eye or with small pupils.
<i>Composition</i>	Trypan blue sterile ophthalmic solution at 0.06% W/V.
<i>Indications</i>	<ul style="list-style-type: none"> • Capsulorhexis. • Phacoemulsification. • Small incision cataract surgery. • Extracapsular cataract surgery. • Anterior/posterior capsulotomy. • Membrane peeling. • Improved visualization.
<i>Sterilization method</i>	Steam.
<i>Supply</i>	Box: 10 units. Contains a 1ml prefilled glass syringe and Luer lock 27G cannula .
<i>Shelf-life</i>	3 years

CE 0434

REV: 02