



The future of vision

ALOS EYELID IMPLANTS



ALOS Eyelid Implants

<i>Brand name</i>	ALOS EYELID IMPLANTS IN GOLD AND PLATINUM						
<i>Indication</i>	Eyelid implants are designed for lagophthalmos functional defect treatment.						
<i>Characteristics</i>	<ul style="list-style-type: none"> • Highly biocompatible materials. • Due to their fine, spherical design, the implants adapt to the curvature of the eyeball, to achieve a better aesthetic result. • Rounded edges to avoid the extrusion of the implant in the long-term. • Suture holes with fixation channels. • Purity of 99.9%. • 7 standard sizes. • Different sizes and customized implants available upon request. • GOLD: Sterile / PLATINUM: Non-sterile. 						
GOLD							
<i>Reference</i>	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18
<i>Weight (g)</i>	0.6	0.8	1.0	1.2	1.4	1.6	1.8
<i>Length (mm)</i>	10.5	12.5	15.5	17.5	19.5	21.5	24.5
<i>Width (mm)</i>	5.7	5.7	5.7	5.7	5.7	5.7	5.7
<i>Thickness (mm)</i>	0.6	0.6	0.6	0.6	0.6	0.6	0.6
<i>Sterilization method</i>	Steam.						
<i>Shelf-life</i>	3 years						

CE 0318

PLATINUM							
<i>Reference</i>	PTH 06	PTH 08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18
<i>Weight (g)</i>	0.6	0.8	1.0	1.2	1.4	1.6	1.8
<i>Length (mm)</i>	10.2	12.4	15.0	17.0	18.8	21.0	24.0
<i>Width (mm)</i>	5.7	5.7	5.7	5.7	5.7	5.7	5.7
<i>Thickness (mm)</i>	0.6	0.6	0.6	0.6	0.6	0.6	0.6
<i>Sterilization method</i>	Non-sterile. Steam sterilization at 121°C for 15 minutes is recommended.						

CE 0434

REV: 02