



The future of vision

EYE SHIELD



Eye Shield

Brand name	AJL EYE SHIELD		
Reference	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
Size	Small: 22 x 19 x 1.7 mm	Medium: 24 x 21 x 1.7 mm	Large: 25 x 22 x 1.7 mm
Product description	AJL EYE SHIELD is a medical device specially designed to protect the eyeball.		
Material / composition	Acrylic.		
Indications	AJL EYE SHIELD is designed to protect the eyeball during eyelid surgery.		
Sterilization method	Ethylene oxide (ETO).		
Supply	Sterile, single container. Device which includes a bulb syringe that facilitates the positioning and extraction of the ocular protector.		
Shelf-life	4.5 years		

CE 0434

REV: 02