

Optec® 5000 Series Vision Tester

Optec® 5000/5000P - Manual Control

Optec® 5500/5500P - Remote Control



5000/5000P Manual Control



5500/5500P Remote Control

Increase profitability by reducing chair time with simple test sequences that can be administered by ancillary office staff in less than 5 minutes

Features include:

- **New, unsurpassed homogenous illumination**
- **L.E.D. lighting system**
- **Background luminance complies with international standards**
- **Distance/Near lens systems**
Two separate viewing ports for accurate distance 20 ft. (6 m) & near 16 in. (40 cm) testing
- Select a standard test package or **customize** to specific testing needs by choosing up to 12 slides from the library of over **150** testing slides
- Standard test packages have been designed for:

Ophthalmology	Public Health
Optometry	Primary Care
Occupational Medicine	D.O.T.
Schools	Driver Rehabilitation
Pediatrics	Contrast Sensitivity (optional)
Athletics	or Customize a Test Package
F.A.A.	
- Optional slide packages include:
 - F.A.C.T.® Contrast Sensitivity (Sine Wave)**
 - C.A.T.® Contrast Sensitivity (Letter)**
- All test slides are produced on high-quality photographic film sealed between 2 glass plates, ensuring that the high-quality images will last for years
- Lightweight and portable, weighing less than 15 lbs (6.8 kg), requiring only 2 square feet (0.6 m) of space, has scratch resistant and washable surface made of ABS plastic

- The Optec® 5500/5500P with remote control includes: accessory kit*, and 1 standard slide package

- The 5500P includes the peripheral test

*Accessory kit includes: dust cover, pointer, lens cleaner, headrest tissues, recording forms, and instruction manual.

- The Optec® 5000/5000P with manual control includes: accessory kit*, and 1 standard slide package

- The 5000P includes the peripheral test