




The future of vision

# PERFORATED CONFORMER



## Perforated Conformer

<i>Brand name</i>	<b>AJL PERFORATED CONFORMER</b>		
<i>Reference</i>	<b>CF-S-AJL</b>	<b>CF-M-AJL</b>	<b>CF-L-AJL</b>
<i>Size</i>	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
<i>Product description</i>	The AJL PERFORATED CONFORMER is a medical device designed to avoid adherence between the tarsal and bulbar conjunctiva. It is transparent and has two orifices for ventilation at the centre.		
<i>Material/ composition</i>	High resistance acrylic material suitable for ocular use.		
<i>Indications</i>	AJL PERFORATED CONFORMER is used to achieve the separation between tarsal and bulbar conjunctiva. The device prevents the shrinkage of the cavity and maintains the orbital space properly after enucleation/visceration.		
<i>Sterilization method</i>	Ethylene Oxide (ETO).		
<i>Supply</i>	Sterile, single container.		
<i>Shelf-life</i>	4.5 years		

CE 0434

REV: 02