



The future of vision

SIOBAL








SIOBAL Silicone Oil

Reference	S1000	S2000	S5000
Brand name	SIOBAL		
Definition	Tamponade agent for retinal detachment.		
Material	High-purity silicone oil free of low-molecular-weight components.		
Viscosity	1 000 cSt	2 000 cSt	5 000 cSt
Indications	<ul style="list-style-type: none"> • Retinal detachment. • Proliferative vitreoretinopathy. • Perforating wounds. • Proliferative diabetic retinopathy. 		
Sterilization method	Steam.		
Supply	Sterile, single unit. 10ml prefilled glass syringe.		
Shelf-life	3 years		

CE 0318

REV: 03

SIOBAL PACK






Reference	PS1000-A* PS1000-OA* PS1000-OR* PS1000-B*		PS2000-A* PS2000-OA* PS2000-OR* PS2000-B*		PS5000-A* PS5000-OA* PS5000-OR* PS5000-B*	
Brand name	SIOBAL PACK: SIOBAL Silicone Oil + Injection Tubing Set (Sterile)					
Viscosity- oil	1000 cSt		2000 cSt		5000 cSt	
Material - oil	High-purity silicone oil free of low-molecular-weight components					
Content	Cannula 20G	Cannula 23G	Cannula 25G	Connection system for injection / extraction	10ml oil prefilled syringe	
						
Supply	5 packs per box. Products sterilized and packed in their original containers by the manufacturer.					
Quality assurance	CE 1014			CE 0318		

***Selection of compatible connection system for the following equipments:**

(A) Alcon: Accurus, Constellation | Optikon: (OA) Assistant, (OR) R-Evolution. | (B) Bausch&Lomb: Millenium, Stellaris



SIOBAL SET

Refnce	ST-A*		ST-OA*		ST-OR*		ST-B*	
Brand name	SIOBAL SET: Silicone Oil Injection / Extraction Tubing Set (Sterile)							
Content	Cannula 20G	Cannula 23G	Cannula 25G	Connection system for injection / extraction	Syringe 10ml			
								
Supply	5 sets per box. Products sterilized and packed in their original containers by the manufacturer.							
Quality assurance	CE 1014			CE 0318		CE 0050		

***Selection of compatible connection system for the following equipments:**

(A) Alcon: Accurus, Constellation | Optikon: (OA) Assistant, (OR) R-Evolution. | (B) Bausch&Lomb: Millenium, Stellaris

REV: 03