



The future of vision

SYMBLEPHARON RING



Syம்பlepharon Ring

<i>Brand name</i>	AJL SYMBLEPHARON RING		
<i>Reference</i>	SF-S-AJL	SF-M-AJL	SF-L-AJL
<i>Size</i>	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
<i>Product description</i>	AJL SYMBLEPHARON RING is a medical device designed with a geometry specifically adapted to the eyeball.		
<i>Material</i>	Acrylic.		
<i>Indications</i>	AJL SYMBLEPHARON RING is used to avoid total or partial adhesion between the internal face of the eyelid and the eyeball, as well as maintaining the depth of the conjunctival bag.		
<i>Sterilization method</i>	Ethylene Oxide (ETO).		
<i>Supply</i>	Sterile, single container.		
<i>Shelf-life</i>	4.5 years		

CE 0434

REV: 02