

Combi 100

Compact refraction unit

Right version

REF 8105-0101-01

Left version

REF 8105-0102-01

If limitations in space for a refraction unit make the use of a small refraction unit essential Combi 100 will be applied.

Features

In order to conduct a complete and of course competent refraction, Combi 100 offers a two-instrument table as well as the fittings to mount a phoropter arm to hold a manual or automatic phoropter (optional). The projector is mounted on the aluminium column.

Small, functional and convenient

- Two-instrument table
- Integrated control panel.

Standard refraction chair

- Stepp-less electrical high adjustment
- Footrest.

Technical specifications

General	
Min. seat height	530 mm
Max. seat height	720 mm
Inclination	Not tiltable
Seat adjustment	No
Arm rest	Can be swivelled
Standards	MDD, CE

