

**A.R.C.
LASER**

enlighten your surgery.

WOLF
SURGICAL LASER


TruBlue



**PHOTO-
ANGIOLYTIC**

CUTS LIKE CO₂

**FIBER GUIDED
LASER**



WOLF TruBlue Laser


ADVANTAGES:

- best coagulation for bloodless surgery
- non-contact cutting, low necrosis zone
- minimal bending radius to target any treatment area
- widest range of application in a „shoe box“





The smallest “**CO₂ + KTP**”-like laser unit for coagulation, vaporisation and bloodless surgery, with all the benefits of a razor sharp scalpel inside an extremely flexible glass fiber. Join the TruBlue revolution!




Versatility meets Efficiency

The Wolf TruBlue gives you a vast versatility in applications which no laser has reached so far. You are also free to apply the laser in contact or non-contact mode to meet the most beneficial requirements of your desired application in

- Phonosurgery
- Laryngology
- Otology

Truly Convincing: TruBlue allows to cut like a CO₂ laser whilst coagulating like a KTP laser – thanks to its peak absorption in hemoglobin and its permeability in water.



TruBlue – Don't call it CO₂

The shockproof and portable WOLF offers unprecedented design, size and function. The laser is supporting a variety of fibers with 300, 400 or 600 micron, providing you with high flexibility. Advances in fiber optics have made them ideally suited for endoscopic use, to travel through routes in the human body where no hand or scalpel has gone before.

Advantages

- Convincing wavelength properties and flexibility.
- Allows outpatient treatment – anywhere.
- Lower costs of accessories and maintenance.

Superb Quality and Safety

For us quality and safety come first. The unit offers superb quality and is backed by a full 2 year warranty.

