


Silkann[®]

PERFECTION

for injections




Sterimedix[®]
aesthetic

Silkann[®]

Aesthetic injections
need a cannula that is reliable,
consistent and cost effective:
perfect in every way.


Silkann cannulas are made to the highest standards.



Silkann cannulas feature a unique polycarbonate hub, with a screw thread.




- Indicator on the hub to show the port orientation
- High quality packaging
- Manufactured in the UK by Sterimedix.

All Silkann products are conveniently packed in a blister containing both the cannula and pre-hole needle





DERMAL FILLER CANNULAS

22 GAUGE CANNULA + 21 GAUGE PRE-HOLE NEEDLE


			
AJ 1831A	22g TW x 70mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 21g pre-hole needle	Box 25
			
AJ 1832A	22g TW x 60mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 21g pre-hole needle	Box 25
			
AJ 1833A	22g TW x 50mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 21g pre-hole needle	Box 25

23 GAUGE CANNULA + 21 GAUGE PRE-HOLE NEEDLE

			
AJ 1840A	23g UTW x 70mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 21g pre-hole needle	Box 25
			
AJ 1842A	23g UTW x 50mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 21g pre-hole needle	Box 25

NEW GTI CANNULA®


25 GAUGE GTI CANNULA + 23 GAUGE PRE-HOLE NEEDLE

			
AJ 1852A	25g UTW x 40mm	GTI Cannula® (Amar) + 23g pre-hole needle	Box 10

Turn to back page for more information about the new GTI Cannula®.

FAT TRANSFER CANNULAS



25 GAUGE FAT TRANSFER CANNULA

			
AE 1702	18g UTW x 70mm	Straight fat transfer cannula with screw thread hub, closed end and single side port	Box 40


20 GAUGE FAT TRANSFER CANNULA

			
AE 1711	20g UTW x 60mm	Straight fat transfer cannula with screw thread hub, closed end and single side port	Box 40


25 GAUGE CANNULA + 23 GAUGE PRE-HOLE NEEDLE

			
AJ 1850A	25g UTW x 50mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 23g pre-hole needle	Box 25
			
AJ 1851A	25g UTW x 40mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 23g pre-hole needle	Box 40


27 GAUGE CANNULA + 25 GAUGE PRE-HOLE NEEDLE

			
AJ 1861A	27g UTW x 40mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 25g pre-hole needle	Box 40

30 GAUGE CANNULA + 27 GAUGE PRE-HOLE NEEDLE





			
AJ 1870A	30g UTW x 30mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port + 27g pre-hole needle	Box 40

18 GAUGE CANNULA

			
AJ 1701	18g UTW x 70mm	Straight dermal cannula with fine thread, two-start hub, closed end and single side port	Box 40


SHARP NEEDLES

SHARP NEEDLES

			
M0286	27g TW x 12.7mm	Straight dermal injection needle	Box 100
			
M0288	29g TW x 12.7mm	Straight dermal injection needle	Box 100
			
M0287	30g TW x 12.7mm	Straight dermal injection needle	Box 100
			
M0289A	30g TW x 4mm	Straight dermal injection needle	Box 100

The latest innovation for small scars and surface defects.

The new Silkann® GTI Cannula®



- 1: The Grooved Tipped Injection Cannula is designed to move through higher density & fibrotic tissue in the face.
- 2: Micro subcision of fibrous strands.
- 3: Fibrous strands are released.
- 4: Injection of dermal filler to fill the void upon cannula withdrawal.

Designed in conjunction with DR Olivier Amar; an Aesthetic & Reconstructive Surgeon based in London, UK.



NEW GTI CANNULA®

25 GAUGE GTI CANNULA + 23 GAUGE PRE-HOLE NEEDLE



AJ 1852A	25g UTW x 40mm	GTI Cannula® (Amar) + 23g pre-hole needle	Box 10
----------	----------------	---	-----------

Distributed by:

For more information call:
+44 (0)1527 501480
or email: sales@sterimedix.com

Sterimedix® aesthetic

Sterimedix Limited, Thornhill Road, North Moons Moat
Redditch, Worcestershire B98 9ND, UK

sterimedix.com