


On Request

## Donnenfeld


Micro soft IOL cutter  
4 mm long blades  
micro notches in blades for traction  
**AA 3749 / AAT 3749**



On Request


## Donnenfeld

Micro IOL holding forceps  
21 G blunt, paddle-shaped tips  
central hole for traction and visibility  
for grasping a soft IOL during cutting and explanation.  
**AA 3749A / AAT 3749A**




## Kreshner

Micro incision  
capsulorhexis forceps, curved  
**AA 3871 / AAT 3871** 20 G  
**AA 3872 / AAT 3872** 23 G



## Vitreous Forceps

Smooth jaws, straight  
**AA 3960 / AAT 3960** 20 G  
**AA 3960A / AAT 3960A** 23 G




## Diamondized Gripping

Forceps straight

**AA 3961 / AAT 3961** 20 G


**AA 3961A / AAT 3961A** 23 G

**AA 3961C / AAT 3961C** 27 G



## Vitreous Forceps

**AA 3962 / AAT 3962** 18 G



## Crocodile Gripping

Forceps

**AA 3963 / AAT 3963** 20 G

**AA 3963A / AAT 3963A** 23 G




## Vitreous Forceps

Serrated jaws, straight

**AA 3965 / AAT 3965** 20 G


**AA 3965A / AAT 3965A** 23 G



## Vitreous Forceps

Smooth jaws, 1x2 teeth straight


**AA 3970 / AAT 3970** 20 G



## End Gripping Forceps

**AA 3971 / AAT 3971** 20 G

**AA 3971A / AAT 3971A** 23 G




## End Gripping Forceps

With micro jaws,

**AA 3972 / AAT 3972** 20 G

**AA 3972A / AAT 3972A** 23 G




## Vitreous Forceps

Smooth jaws, 0.3 mm platform, straight

**AA 3975 / AAT 3975** 20 G


**AA 3975A / AAT 3975A** 23 G



## Pick Forceps

**AA 3976 / AAT 3976** 20 G


**AA 3976A / AAT 3976A** 23 G



## Appa Pick Forceps

**AA 3977 / AAT 3977** 20 G

**AA 3977A / AAT 3977A** 23 G



## Gripping Angled Forceps

**AA 3978 / AAT 3978** 20 G



## Sub Retinal Forceps

**AA 3979 / AAT 3979** 20 G



## Vitreous Forceps

Smooth jaws, vertical

**AA 3980 / AAT 3980** 20 G

**AA 3980A / AAT 3980A** 23 G




## Vitreous Forceps

30° Ang jaws

**AA 3981 / AAT 3981** 23 G

**AA 3981B / AAT 3981B** 27 G



## Asymmetrical

End gripping forceps

AA 3982 / AAT 3982 20 G

AA 3982A / AAT 3982A 23 G



## Tanaka

MaCulorhexis forceps

AA 3983 / AAT 3983 23 G



## Vitreous

Forceps


AA 3985 / AAT 3985 20 G



## ILM Forceps

AA 3986 / AAT 3986 20 G

AA 3986A / AAT 3986A 23 G




## Lens Cartridge

Loading forceps ICL

AA 3987 / AAT 3987 20 G

AA 3987A / AAT 3987A 23 G

AA 3987C / AAT 3987C 27 G



## Foreign Body Forceps

AA 3988 / AAT 3988 17 G



## Vitreoretinal Forceps

With cup jaws

AA 3989 / AAT 3989 20 G



## Iridectomy Forceps

Designed to perform capsulorhexis through 0.8 mm side-port incision


straight 23 gauge shaft

delicate tips

with lock

compatible with squeeze handle


AA 3990 / AAT 3990



## Spring Gripping Forceps


AA 3991 / AAT 3991 20 G

AA 3991A / AAT 3991A 23 G




## Large Foreign Body Forceps

AA 3992 / AAT 3992 20 G



## Endothelial Gripping Forceps

AA 3993 / AAT 3993 20 G



## Agarwal Glued IOL Forceps

AA 3995 / AAT 3995 straight 20 G

AA 3995A / AAT 3995A straight 23 G

AA 3996 / AAT 3996 curved 20 G

AA 3996A / AAT 3996A curved 23 G



On Request

## Appasamy

Intraocular scissors  
squeeze handle, protective cap  
17gauge to cut capsular tags

**AA 4400** straight

**AA 4410** curved

blade length 7 mm



On Request

## Scimeca

Lens loop amputator  
for cutting poly propylene  
or PMMA lens haptics in the  
anterior or posterior chamber

**AA 4420** left

**AA 4425** right

blade length 3 mm



## Vitreous

Scissors  
straight jaws

blade length 5 mm

**AA 4450 / AAT 4450** 20 G

blade length 4 mm

**AA 4450A / AAT 4450A** 23 G



## Vitreous

Scissors  
Curved jaws

blade length 5 mm

**AA 4455 / AAT 4455** 20 G

blade length 4 mm

**AA 4455A / AAT 4455A** 23 G



## Vitreous

Scissors

angled jaws 30°

blade length 3 mm

**AA 4460 / AAT 4460** 20 G



## Vitreous


Scissors

vertical cut 60°

blade length 3 mm

**AA 4470 / AAT 4470** 20 G

**AA 4470A / AAT 4470A** 23 G



## Vitreous

Scissors

vertical cut 90°

blade length 3 mm

**AA 4480 / AAT 4480** 20 G

**AA 4480A / AAT 4480A** 23 G



## Vertical Scissors

55° 20 G

**AA 4605 / AAT 4605**



## Horizontal Scissors

45°, regular blades

**AA 4610 / AAT 4610** 20 G



## Horizontal Scissors

45°, short blades

**AA 4615 / AAT 4615** 20 G



## Horizontal Scissors

Sub Retinal scissors, curved

**AA 4620 / AAT 4620** 23 G

