


The latest innovation for small scars and surface defects.

The new Silkann® GTI Cannula®

Designed specially to disrupt higher density fibrous scar tissue to perform a simple and safe subcision.




- 1: The Grooved Tipped Injection Cannula is designed to move safely through higher density & fibrotic tissue in the face.
- 2: Micro subcision of fibrous strands.
- 3: Fibrous strands are released.
- 4: Injection of dermal filler to fill the void upon cannula withdrawal.

Designed in conjunction with DR Olivier Amar; an Aesthetic & Reconstructive Surgeon based in London, UK.



NEW GTI CANNULA®

25 GAUGE GTI CANNULA + 23 GAUGE PRE-HOLE NEEDLE



AJ 1852A	25g UTW x 40mm	GTI Cannula® (Amar) + 23g pre-hole needle	Box 10
----------	----------------	---	-----------

Distributed by:

For more information call:
+44 (0)1527 501480
or email: sales@sterimedix.com

Sterimedix®
aesthetic

Sterimedix Limited, Thornhill Road, North Moons Moat
Redditch, Worcestershire B98 9ND, UK

sterimedix.com