



ESNOOPER



The future of vision

# ESNOOPER Glaucoma Implant

Reference	V2000	Brand name	ESNOOPER
Definition	ESNOOPER V2000 Scleral Implant.		
Indication	The ESNOOPER V2000 scleral implant is indicated for non-penetrating deep sclerectomy (NPDS) in open-angle glaucoma surgery.		
Product description	The ESNOOPER V2000 scleral implant is a non-reabsorbable implant that provides a fixed, permanent intrascleral space that facilitates the drainage of aqueous humour.		
Dimensions	3.00mm X 1.40mm X 2.85mm.		
Orifices	2 (suture/drainage).		
Longitudinal orifice	Drainage.		
Material	Hydroxyethyl Methacrylate (HEMA).		
Sterilisation method	Steam.		
Supply	The ESNOOPER V2000 implant comes in an individual sterile vial housed in a clamp support and immersed in a liquid medium to maintain proper device hydration.		
Quality Assurance	CE 0434 marking.	Shelf-life	3 years

## EPNP surgical equipment

Definition	EPNP Kit (non-penetrating deep sclerectomy).		
Brand name	KIT EPNP	Product code	AJL- KIT-001
Description	<p>The EPNP Kit is presented in a rigid plastic blister in which the following references are grouped (*):</p> <ul style="list-style-type: none"> <li>• <b>EPNP-2001</b> Liebermann Blepharo 14mm.</li> <li>• <b>EPNP-2002</b> Double Scleral Marker.</li> <li>• <b>EPNP-2003</b> Mermoud Spatula.</li> <li>• <b>EPNP-2004</b> Mermoud Clamp.</li> <li>• <b>EPNP-2005</b> Curved micro-scissor.</li> <li>• <b>PST-4600A</b> Sterilisation box (10x16x2 cm).</li> </ul>		
Sterilisation method	Non-sterile.	Quality Assurance	CE marking.



\* The instruments can be supplied as a set (EPNP Kit) or separately.

REV: 02