

**REFRESH &
REJUVENATE**
DRY EYES
NATURALLY

Bruder
EyeCompress
MOIST HEAT



The problem

Dry, irritated, itchy eyes

Let your doctor know if you have any of the following symptoms:

- Dry, itchy or burning eyes
- Watery, gritty or blurry eyes
- Sensitivity to light, causing you to squint or blink
- Difficulty focusing due to dryness

Dryness is caused by a variety of factors:

- Age
- Eye surgery
- Extended computer use
- Contact lens wear
- Certain medications
- Allergies
- Eye disorders
- Diseases (Arthritis, Diabetes, etc.)
- Hormonal imbalance
- Environmental factors

Contact us:

(888) 827-8337

A simple & natural solution

The **Bruder Moist Heat Eye Compress** delivers a doctor-recommended moist heat treatment in just seconds from the microwave.



*Bruder Moist Heat Eye Compress with MediBeads
ITEM #34160*

Ask your doctor why adding the easy-to-use **Bruder Moist Heat Eye Compress** to your eye wellness routine will help relieve your symptoms safely, effectively, and economically.

The solution

Bruder Moist Heat Eye Compress

The compress relieves dryness simply, naturally, and conveniently. Just place the compress in a clean microwave for 20-30 seconds, then apply over closed eyes for 10 minutes. The **moist heat** provides welcome relief in minutes. The compress will maintain prescribed heat levels for the entire treatment. Use as often as necessary, typically one or two times a day.

The **Bruder Moist Heat Eye Compress** will not interfere with other prescriptions or over-the-counter medications. It can be hand washed when needed and will last 6-9 months with frequent daily use.

How it works:

The **Bruder Moist Heat Eye Compress** features patented MediBeads® technology. Water molecules are continuously absorbed from the air and stored inside the beads. When microwaved, the absorbed water is released as clean moist heat. There is no wait, no mess, and no need to add water.

The moist heat provided by the Bruder Moist Heat Eye Compress may also be a perfect solution for treating Blepharitis, Styes, and Chalazions.

Available from:



**Ask your doctor how Bruder products
can benefit you.**

#34160 — Bruder Moist Heat Eye Compress

#34170 — Bruder Moist Heat Eye Compress - Single Eye

#34180 — Bruder Moist Heat Sinus Compress

#61080 — Bruder Cold Therapy Eye Compress

#78900 — Bruder Hygienic Eyelid Sheets

#46201 — Bruder Hygienic Eyelid Solution

Bruder[®]

© Bruder Healthcare Company, LLC

1115 Ridgeland Parkway, #101 Alpharetta, GA 30004

Bruder and MediBeads are trademarks of
Bruder Healthcare Company.

To learn more:

(888) 827-8337

bruder.com

eyes@bruder.com

#B34160P_V1903